

Anti-Pax9 Monoclonal Antibody For Wb, Elisa - P55771

Item Number: CM0003172

Introduction

Anti-PAX9 rabbit monoclonal antibody validated for WB and ELISA. Detects human, mouse, rat PAX9 transcription factor. Ideal for developmental biology and transcription research.

[Learn More](#)

Parameter	Value
Product Name	PAX9 Rabbit mAb
Remarks/Alias	STHAG3; PAX9
Species	Human
Gene ID (Human)	5083
Gene ID	5083
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PAX9 (P55771)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55771
Protein Weight	36kDa
Shipping	Ice bag