


Anti-Paxillin Monoclonal Antibody For Wb, Elisa - P49023

Item Number: CM0003674



Introduction

Recombinant rabbit monoclonal antibody against human Paxillin, validated in WB and ELISA. Reacts with human, mouse, and rat samples. Investigates focal adhesion, cytoskeleton organization, and cell junction dynamics.

[Learn More](#)

Parameter	Value
Product Name	Paxillin Rabbit mAb
Remarks/Alias	Paxillin; PXN; paxillin
Species	Human
Gene ID (Human)	5829
Gene ID	5829
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Paxillin (P49023).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49023
Protein Weight	65kDa/61kDa/66kDa/47kDa
Shipping	Ice bag