


## Anti-Paxillin Rabbit Polyclonal Antibody For Wb, Elisa - P49023

Item Number: CM0021095



### Introduction

Rabbit polyclonal antibody targeting human Paxillin, suitable for Western blot and ELISA. Detects Paxillin in human and mouse samples for studying focal adhesion dynamics and cytoskeletal interactions.

[Learn More](#)

| Parameter        | Value                                                                                                  |
|------------------|--------------------------------------------------------------------------------------------------------|
| Product Name     | Paxillin Rabbit pAb                                                                                    |
| Remarks/Alias    | Paxillin; PXN; paxillin                                                                                |
| Species          | Human                                                                                                  |
| GeneID (Human)   | 5829                                                                                                   |
| GeneID           | 5829                                                                                                   |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Paxillin (NP_002850.2) |
| Source           | Rabbit                                                                                                 |
| Category         | Polyclonal Antibodies                                                                                  |
| Application      | WB, ELISA                                                                                              |
| Cross Reactivity | Human, Mouse                                                                                           |
| SWISS            | P49023                                                                                                 |
| Protein Weight   | 65kDa                                                                                                  |
| Shipping         | Ice bag                                                                                                |