

Anti-Pbef/Visfatin/Nampt Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P43490

Item Number: CM0016391



Introduction

Rabbit polyclonal antibody targeting human PBEF/Visfatin/NAMPT. Validated for WB, IHC-P, IF/ICC, and ELISA with cross reactivity to human, mouse, and rat. Suitable for NAD metabolism, circadian rhythm, and cytokine studies.

[Learn More](#)

Parameter	Value
Product Name	PBEF/Visfatin/NAMPT Rabbit pAb
Remarks/Alias	VF; PBEF; PBEF1; VISFATIN; 1110035O14Rik; PBEF/Visfatin/NAMPT
Species	Human
GeneID (Human)	10135
GeneID	10135
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43490
Protein Weight	56kDa
Shipping	Ice bag