

Anti-Pbef/Visfatin/Nampt Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P43490

Item Number: CM0007190



Introduction

High-performance rabbit monoclonal antibody against PBEF/Visfatin/NAMPT (UniProt P43490). Validated for Western blot, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat samples. Ideal for NAD biosynthesis, immunomodulation, and adipokine research.

[Learn More](#)

Parameter	Value
Product Name	PBEF/Visfatin/NAMPT Rabbit mAb
Remarks/Alias	VF; PBEF; PBEF1; VISFATIN; 1110035O14Rik; PBEF/Visfatin/NAMPT
Species	Human
GeneID (Human)	10135
GeneID	10135
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 300-400 of human PBEF/Visfatin/NAMPT (NP_005737.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43490
Protein Weight	56kDa
Shipping	Ice bag