

## Anti-Pbk Polyclonal Antibody For Wb, Elisa - Q96Kb5

Item Number: CM0021889



### Introduction

Rabbit polyclonal antibody targeting human PBK (TOPK) kinase, validated for WB and ELISA applications. Cross-reacts with mouse and rat. Suitable for cell cycle, DNA damage response, and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                          |
|------------------|----------------------------------------------------------------------------------------------------------------|
| Product Name     | [KO Validated] PBK Rabbit pAb                                                                                  |
| Remarks/Alias    | SPK; CT84; TOPK; HEL164; Nori-3; PBK                                                                           |
| Species          | Human                                                                                                          |
| Gene ID (Human)  | 55872                                                                                                          |
| Gene ID          | 55872                                                                                                          |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-322 of human PBK (NP_060962.2) |
| Source           | Rabbit                                                                                                         |
| Category         | Polyclonal Antibodies                                                                                          |
| Application      | WB, ELISA                                                                                                      |
| Cross Reactivity | Human, Mouse, Rat                                                                                              |
| SWISS            | Q96KB5                                                                                                         |
| Protein Weight   | 36kDa                                                                                                          |
| Shipping         | Ice bag                                                                                                        |