


Anti-Pbrm1 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q86U86

Item Number: CM0002285



Introduction

Rabbit monoclonal antibody targeting PBRM1 (BAF180), validated for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, rat. Ideal for chromatin remodeling and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	PBRM1 Rabbit mAb
Remarks/Alias	PB1; RCC; BAF180; SMARCH1; PBRM1
Species	Human
Gene ID (Human)	55193
Gene ID	55193
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 530-680 of human PBRM1 (Q86U86)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86U86
Protein Weight	193 kDa
Shipping	Ice bag