


Anti-Pbrm1 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q86U86

Item Number: CM0027261



Introduction

Rabbit polyclonal antibody against human PBRM1 (BAF180), validated for WB, IP, and ELISA. Cross-reacts with mouse and rat. Targets a SWI/SNF chromatin remodeling factor with tumor suppressor activity.

[Learn More](#)

Parameter	Value
Product Name	PBRM1 Rabbit pAb
Remarks/Alias	PB1; RCC; BAF180; SMARCH1; PBRM1
Species	Human
Gene ID (Human)	55193
Gene ID	55193
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human PBRM1 (NP_060783.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86U86
Protein Weight	193kDa
Shipping	Ice bag