

Anti-Pbr/Tspo Monoclonal Antibody For Wb - P30536

Item Number: CM0009004

Introduction

Recombinant rabbit monoclonal antibody against human PBR/TSPO (TSPO, PBR), validated for WB, IHC-P, IF/ICC, ELISA. Detects human and mouse protein. Ideal for mitochondrial and cholesterol transport research.

[Learn More](#)

Parameter	Value
Product Name	PBR/TSPO Rabbit mAb
Remarks/Alias	DBI; IBP; MBR; PBR; PBS; BPBS; BZRP; PKBS; PTBR; mDRC; pk18; TSPO1; PBR/TSPO
Species	Human
GeneID (Human)	706
GeneID	706
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 70-169 of human PBR/TSPO (P30536).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P30536
Protein Weight	19kDa
Shipping	Ice bag