

Anti-Pbx1 Polyclonal Antibody For Wb, Elisa - P40424

Item Number: CM0028712



Introduction

High-quality rabbit polyclonal antibody against human PBX1 transcription factor. Validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Suitable for research in transcriptional regulation, development, and cancer.

[Learn More](#)

Parameter	Value
Product Name	PBX1 Rabbit pAb
Remarks/Alias	CAKUHD; PBX1
Species	Human
GeneID (Human)	5087
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 311-430 of human PBX1 (NP_002576.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40424
Protein Weight	47kDa
Shipping	Ice bag