

Anti-Pbx3 Polyclonal Antibody For Wb And Elisa - P40426

Item Number: CM0024601

Introduction

PBX3 Rabbit pAb is a polyclonal antibody recognizing human PBX3 transcription factor. Validated for WB and ELISA, it cross-reacts with human and mouse. Ideal for transcription regulation and developmental studies.

[Learn More](#)

Parameter	Value
Product Name	PBX3 Rabbit pAb
Remarks/Alias	PBX3
Species	Human
Gene ID (Human)	5090
Gene ID	5090
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PBX3 (NP_006186.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P40426
Protein Weight	47kDa
Shipping	Ice bag