

Anti-Pcbd1 Monoclonal Antibody For Wb, Ihc-P, Elisa - P61457

Item Number: CM0009664

Introduction

High-quality rabbit monoclonal antibody against human PCBD1 (Pterin-4-alpha-carbinolamine dehydratase). Validated for WB, IHC-P, ELISA; reacts with human, mouse, rat. Suitable for tetrahydrobiopterin biosynthesis and transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	PCBD1 Rabbit mAb
Remarks/Alias	PCD; PHS; DCOH; PCBD; PCBD1
Species	Human
GeneID (Human)	5092
GeneID	5092
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-104 of human PCBD1 (P61457).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61457
Protein Weight	12kDa
Shipping	Ice bag