

Anti-Pcbd1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P61457

Item Number: CM0028321

Introduction

Rabbit polyclonal antibody against human PCBD1, validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for tetrahydrobiopterin biosynthesis and transcription research.

[Learn More](#)

Parameter	Value
Product Name	PCBD1 Rabbit pAb
Remarks/Alias	PCD; PHS; DCOH; PCBD; PCBD1
Species	Human
Gene ID (Human)	5092
Gene ID	5092
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61457
Protein Weight	12kDa
Shipping	Ice bag