


Anti-Pcca Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P05165

Item Number: CM0029041



Introduction

Rabbit polyclonal anti-PCCA antibody validated for western blot, immunohistochemistry-paraffin, and ELISA. Cross-reacts with human, mouse, and rat. Targets propionyl-CoA carboxylase alpha chain, involved in branched-chain amino acid catabolism.

[Learn More](#)

Parameter	Value
Product Name	PCCA Rabbit pAb
Remarks/Alias	PCCA
Species	Human
GeneID (Human)	5095
GeneID	5095
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 449-728 of human PCCA (NP_000273.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05165
Protein Weight	80kDa
Shipping	Ice bag