


Anti-Pccb Rabbit Monoclonal Antibody [Ko Validated] For Wb, Ihc-P, If/Icc, Elisa - P05166

Item Number: CM0006665



Introduction

KO-validated anti-PCCB rabbit monoclonal antibody for Western blot, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat PCCB, a mitochondrial enzyme involved in propionyl-CoA carboxylation. Research-use-only IVD raw material.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PCCB Rabbit mAb
Remarks/Alias	PCCB; CB
Species	Human
GeneID (human)	5096
GeneID	5096
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 260-539 of Human PCCB [NP_000523.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05166
Protein Weight	58kDa
Shipping	Ice bag