


Anti-Pccb Rabbit Polyclonal Antibody For Wb And Elisa - P05166

Item Number: CM0025095



Introduction

High-quality rabbit polyclonal antibody targeting human propionyl-CoA carboxylase beta chain (PCCB). Validated for Western blotting (WB) and ELISA applications, with cross-reactivity to mouse. Ideal for mitochondrial metabolism, branched-chain amino acid catabolism, and propionic acidemia research.

[Learn More](#)

Parameter	Value
Product Name	PCCB Rabbit pAb
Remarks/Alias	PCCB
Species	Human
GeneID (Human)	5096
GeneID	5096
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 260-539 of human PCCB (NP_000523.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P05166
Protein Weight	58kDa
Shipping	Ice bag