

Anti-Pcdh1 Polyclonal Antibody For Wb And Elisa - Q08174

Item Number: CM0017294



Introduction

Rabbit polyclonal antibody targeting human Protocadherin-1 (PCDH1, PC42) for Western blotting and ELISA. Suitable for cell adhesion and neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	PCDH1 Rabbit pAb
Remarks/Alias	PC42; PCDH42; PCDH1
Species	Human
Gene ID (Human)	5097
Gene ID	5097
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-220 of human PCDH1 (NP_002578.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q08174
Protein Weight	115kDa
Shipping	Ice bag