


Anti-Pcdh15 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96Qu1

Item Number: CM0018499



Introduction

Rabbit polyclonal antibody against human protocadherin-15, validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Targets a protein critical for retinal and cochlear function; suitable for hearing and vision research.

[Learn More](#)

Parameter	Value
Product Name	PCDH15 Rabbit pAb
Remarks/Alias	USH1F; CDHR15; DFNB23; PCDH15
Species	Human
Gene ID (Human)	65217
Gene ID	65217
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 160-400 of human PCDH15 (NP_001136235.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96QU1
Protein Weight	216kDa
Shipping	Ice bag