

Anti-Pcdh9 Polyclonal Antibody For Wb, Elisa - Q9Hc56

Item Number: CM0015593



Introduction

Rabbit polyclonal antibody against human Protocadherin-9 (PCDH9) for Western blot and ELISA applications. Cross-reacts with human samples. Predicted molecular weight 136 kDa.

[Learn More](#)

Parameter	Value
Product Name	PCDH9 Rabbit pAb
Remarks/Alias	PCDH9
Species	Human
GeneID (Human)	5101
GeneID	5101
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 400-600 of human PCDH9 (NP_982354.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9HC56
Protein Weight	136kDa
Shipping	Ice bag