

Anti-Pcdha13 Polyclonal Antibody For Wb, Elisa - Q9Y5I0

Item Number: CM0029068

Introduction

High-quality rabbit polyclonal antibody targeting protocadherin alpha-13 (PCDHA13), validated for Western blot and ELISA. Reacts with human, mouse, and rat. Ideal for neuroscience research on neuronal connections and cell adhesion.

[Learn More](#)

Parameter	Value
Product Name	PCDHA13 Rabbit pAb
Remarks/Alias	CNR5; CNRN5; CNRS5; CRNR5; PCDH-ALPHA13; PCDHA13
Species	Human
Gene ID (Human)	56136
Gene ID	56136
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-360 of human PCDHA13 (NP_061727.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5I0
Protein Weight	102kDa
Shipping	Ice bag