

Anti-Pcdhb16 Polyclonal Antibody For Wb, Elisa - Q9Nrj7

Item Number: CM0022809

Introduction

Rabbit polyclonal antibody against human PCDHB16, validated for WB and ELISA. Cross-reacts with mouse and rat. Useful for studying neuronal cell adhesion and brain-specific connections.

[Learn More](#)

Parameter	Value
Product Name	PCDHB16 Rabbit pAb
Remarks/Alias	ME1; PCDH3X; PCDHB8a; PCDH-BETA16; PCDHB16
Species	Human
Gene ID (Human)	57717
Gene ID	57717
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9NRJ7
Protein Weight	85kDa
Shipping	Ice bag