

Anti-Pcdhb4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y5E5

Item Number: CM0021781

Introduction

Rabbit polyclonal antibody targeting human PCDHB4 (protocadherin beta-4). Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Suitable for cell adhesion and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	PCDHB4 Rabbit pAb
Remarks/Alias	PCDH-BETA4; PCDHB4
Species	Human
Gene ID (Human)	56131
Gene ID	56131
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-300 of human PCDHB4 (NP_061761.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5E5
Protein Weight	87kDa
Shipping	Ice bag