


Anti-Pcf11 Rabbit Polyclonal Antibody For Wb, Elisa - 094913

Item Number: CM0020220



Introduction

Rabbit polyclonal antibody targeting human PCF11, a pre-mRNA cleavage complex II component. Validated for Western blot and ELISA, with cross-reactivity to mouse. Ideal for transcription termination research.

[Learn More](#)

Parameter	Value
Product Name	PCF11 Rabbit pAb
Remarks/Alias	PCF11
Species	Human
Gene ID (Human)	51585
Gene ID	51585
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1256-1555 of human PCF11 (NP_056969.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	094913
Protein Weight	173kDa
Shipping	Ice bag