

Anti-Pcgf1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Bsm1

Item Number: CM0029948



Introduction

Rabbit polyclonal antibody against human PCGF1 (NSPc1, RNF68) validated for WB, IF/ICC, and ELISA. Detects human, mouse, rat. Useful for studying Polycomb group complexes, transcriptional repression, chromatin remodeling, and cell proliferation.

[Learn More](#)

Parameter	Value
Product Name	PCGF1 Rabbit pAb
Remarks/Alias	NSPC1; RNF68; RNF3A-2; 2010002K04Rik; PCGF1
Species	Human
GeneID	84759
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 160-259 of human PCGF1 (NP_116062.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BSM1
Protein Weight	30kDa
Shipping	Ice bag