

Anti-Pcgf6 Rabbit Polyclonal Antibody (Pab) For Wb And Elisa - Q9Bye7

Item Number: CM0028360



Introduction

Rabbit polyclonal antibody against human PCGF6 (MBLR/RNF134), a PRC1 component. Validated for WB and ELISA; cross-reacts with mouse and rat. For epigenetic and developmental studies.

[Learn More](#)

Parameter	Value
Product Name	PCGF6 Rabbit pAb
Remarks/Alias	MBLR; RNF134; PCGF6
Species	Human
Gene ID (Human)	84108
Gene ID	84108
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 161-350 of human PCGF6 (NP_001011663.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYE7
Protein Weight	39kDa
Shipping	Ice bag