

Anti-Pcif1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q9H4Z3

Item Number: CM0006365



Introduction

Rabbit monoclonal antibody against PCIF1, validated for Western blotting, immunohistochemistry (paraffin), immunoprecipitation, and ELISA. Reacts with human, mouse, rat. Suitable for cap methylation studies in mRNA processing.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PCIF1 Rabbit mAb
Remarks/Alias	CAPAM; hCAPAM; hPCIF1; C20orf67; PPP1R121; F1
Species	Human
Gene ID	63935
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human PCIF1(NP_071387.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H4Z3
Protein Weight	81kDa
Shipping	Ice bag