


Anti-Pcif1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9H4Z3

Item Number: CM0018859



Introduction

Rabbit polyclonal antibody targeting human PCIF1 (CAPAM), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Research use for epitranscriptomics and mRNA cap methylation studies.

[Learn More](#)

Parameter	Value
Product Name	PCIF1 Rabbit pAb
Remarks/Alias	CAPAM; hCAPAM; hPCIF1; C20orf67; PPP1R121; PCIF1
Species	Human
GeneID (Human)	63935
GeneID	63935
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human PCIF1 (NP_071387.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H4Z3
Protein Weight	81kDa
Shipping	Ice bag