

Anti-Pcif1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q9H4Z3

Item Number: CM0022997



Introduction

Rabbit polyclonal antibody against human mRNA cap methyltransferase PCIF1 (CAPAM), validated for WB, IHC-P, IP, and ELISA in human, mouse, and rat samples. Targets amino acids 1-200 of NP_071387.1.

[Learn More](#)

Parameter	Value
Product Name	PCIF1 Rabbit pAb
Remarks/Alias	CAPAM; hCAPAM; hPCIF1; C20orf67; PPP1R121; PCIF1
Species	Human
GeneID (Human)	63935
GeneID	63935
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human PCIF1 (NP_071387.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H4Z3
Protein Weight	81kDa
Shipping	Ice bag