

Anti-Pck1 Monoclonal Antibody (Pepck-C) For Western Blot And Elisa - P35558

Item Number: CM0006904



Introduction

Rabbit monoclonal antibody against human PCK1 (PEPCK-C), validated for Western blot and ELISA. Detects human, mouse, and rat PCK1, ideal for gluconeogenesis and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	PCK1 Rabbit mAb
Remarks/Alias	PCKDC; PEPCK1; PEPCKC; PEPCK-C; PCK1
Species	Human
Gene ID	5105
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 523-622 of human PCK1 (NP_002582.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35558
Protein Weight	69kDa
Shipping	Ice bag