

## Anti-Pck1 Rabbit Polyclonal Antibody For Wb - P35558

Item Number: CM0020916



### Introduction

Rabbit polyclonal antibody targeting human PCK1 (PEPCK-C) for WB, IP, IF/ICC, IF-P, and ELISA. Detects endogenous PCK1 in human, mouse, and rat samples. Recombinant immunogen (aa 490-622).

[Learn More](#)

| Parameter        | Value                                                                                                                                  |
|------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PCK1 Rabbit pAb                                                                                                                        |
| Remarks/Alias    | PCKDC; PEPCK1; PEPCKC; PEPCK-C; PCK1                                                                                                   |
| Species          | Human                                                                                                                                  |
| GeneID (Human)   | 5105                                                                                                                                   |
| GeneID           | 5105                                                                                                                                   |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 490-622 of human PCK1 (NP_002582.3). |
| Source           | Rabbit                                                                                                                                 |
| Category         | Polyclonal Antibodies                                                                                                                  |
| Application      | WB, IP, IF/ICC, IF-P, ELISA                                                                                                            |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                      |
| SWISS            | P35558                                                                                                                                 |
| Protein Weight   | 69kDa                                                                                                                                  |
| Shipping         | Ice bag                                                                                                                                |