


Anti-Pcmt1 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P22061

Item Number: CM0028966



Introduction

Rabbit polyclonal antibody against human PCMT1 (PIMT), validated for WB, IF/ICC, IP, and ELISA. Reacts with human, mouse, and rat samples. Ideal for protein repair and aging research.

[Learn More](#)

Parameter	Value
Product Name	PCMT1 Rabbit pAb
Remarks/Alias	PIMT; PCMT1
Species	Human
GeneID (Human)	5110
GeneID	5110
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22061
Protein Weight	25kDa
Shipping	Ice bag