


Anti-Pcna Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P12004

Item Number: CM0003877



Introduction

High-quality anti-PCNA rabbit monoclonal antibody validated for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, rat, and African green monkey. Targets the 29kDa proliferating cell nuclear antigen involved in DNA replication and repair.

[Learn More](#)

Parameter	Value
Product Name	PCNA Rabbit mAb
Remarks/Alias	ATLD2; PCNA
Species	Human
GeneID (Human)	5111
GeneID	5111
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 162-261 of human PCNA (NP_002583.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, African green monkey
SWISS	P12004
Protein Weight	29kDa
Shipping	Ice bag