


Anti-Pcna Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P12004

Item Number: CM0014020



Introduction

Rabbit polyclonal antibody targeting PCNA (P12004), suitable for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for studying DNA replication, repair, and proliferation.

[Learn More](#)

Parameter	Value
Product Name	PCNA Rabbit pAb
Remarks/Alias	ATLD2; PCNA
Species	Human
Gene ID (Human)	5111
Gene ID	5111
Immunogen	Synthetic peptide corresponding to amino acids 162-261 of human PCNA (NP_002583.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12004
Protein Weight	29kDa
Shipping	Ice bag