


Anti-Pcna Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P12004

Item Number: CM0014716



Introduction

Rabbit polyclonal antibody against human PCNA (DNA sliding clamp). Validated for WB, IHC-P, IF/ICC, and ELISA; cross-reacts with mouse and rat. Ideal for DNA replication and repair studies.

[Learn More](#)

Parameter	Value
Product Name	PCNA Rabbit pAb
Remarks/Alias	ATLD2; PCNA
Species	Human
GeneID (Human)	5111
GeneID	5111
Immunogen	Synthetic peptide (aa 162-261 of human PCNA, NP_002583.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12004
Protein Weight	29kDa
Shipping	Ice bag