

Anti-Pcnt Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O95613

Item Number: CM0023653

Introduction

This rabbit polyclonal antibody targeting human PCNT (Pericentrin) is suitable for Western blot, immunofluorescence/immunocytochemistry, and ELISA. It cross-reacts with rat and is an essential tool for centrosome and microtubule research.

[Learn More](#)

Parameter	Value
Product Name	PCNT Rabbit pAb
Remarks/Alias	KEN; PCN; MOPD2; PCNT2; PCNTB; PCTN2; SCKL4; PCNT
Species	Human
Gene ID (Human)	5116
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	O95613
Protein Weight	378kDa
Shipping	Ice bag