

Anti-Pcnx1 Polyclonal Antibody For Wb, Elisa - Q96Rv3

Item Number: CM0020231

Introduction

Rabbit polyclonal antibody targeting human PCNX1 (pecanex-like protein 1). Applicable for Western blot and ELISA. Immunogen: synthetic peptide within amino acids 2250-2341. Validated for human samples.

[Learn More](#)

Parameter	Value
Product Name	PCNX1 Rabbit pAb
Remarks/Alias	PCNX; PCNXL1; pecanex; PCNX1
Species	Human
GeneID (Human)	22990
GeneID	22990
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96RV3
Protein Weight	259kDa
Shipping	Ice bag