


Anti-Pcpe-1 Monoclonal Antibody For Wb And Elisa - Q15113

Item Number: CM0005037



Introduction

High-quality rabbit monoclonal antibody against human PCPE-1, validated for WB and ELISA with reactivity to human, mouse, and rat. Ideal for studies of collagen processing and extracellular matrix biology.

[Learn More](#)

Parameter	Value
Product Name	PCPE-1 Rabbit mAb
Alias	PCPE; PCPE1; PCPE-1
Species	Human
Gene ID	5118
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 304-449 of human TLR4 (NP_002584.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15113
Protein Weight	48 kDa
Shipping	Ice bag