

Anti-Pcsk9 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8Nbp7

Item Number: CM0006871



Introduction

Recombinant anti-human PCSK9 rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects human and mouse PCSK9, essential for cholesterol homeostasis research.

[Learn More](#)

Parameter	Value
Product Name	PCSK9 Rabbit mAb
Remarks/Alias	PCSK9; FH3; HCHOLA3; LDLQC1; NARC-1; NARC1; PC9; proprotein convertase subtilisin/kexin type 9
Species	Human
GeneID (human)	255738
GeneID	255738
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 31-692 of human PCSK9 (NP_777596.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8NBP7
Protein Weight	20 kDa/74 kDa
Shipping	Ice bag