

Anti-Pcsk9 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q8Nbp7

Item Number: CM0021782



Introduction

Rabbit polyclonal antibody against PCSK9, suitable for WB, IF/ICC, IP, and ELISA in human and rat samples. Aiding research in cholesterol homeostasis and LDLR degradation.

[Learn More](#)

Parameter	Value
Product Name	PCSK9 Rabbit pAb
Remarks/Alias	PCSK9; FH3; HCHOLA3; LDLCQ1; NARC-1; NARC1; PC9; proprotein convertase subtilisin/kexin type 9
Species	Human
GeneID	255738
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-692 of human PCSK9 (NP_777596.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8NBP7
Protein Weight	20kDa/74kDa
Shipping	Ice bag