

Anti-Pcsk9 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q80W65

Item Number: CM0019325

Introduction

Rabbit polyclonal antibody against mouse PCSK9 (UniProt Q80W65), validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for cholesterol metabolism and LDL receptor degradation studies.

[Learn More](#)

Parameter	Specification
Product Name	PCSK9 Rabbit pAb
Remarks/Alias	FH3; PC9; Narc1; HCHOLA3; PCSK9
Species	Mouse
GenelD (human)	100102
GenelD	100102
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-694 of mouse PCSK9 (NP_705793.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q80W65
Protein Weight	74KDa
Shipping	Ice bag