

Anti-Pcyt2 Polyclonal Antibody For Wb, Elisa - Q99447

Item Number: CM0015103



Introduction

PCYT2 rabbit polyclonal antibody, validated for Western blot (WB) and ELISA. Detects endogenous PCYT2 in human, mouse, and rat samples. Useful for studying phosphatidylethanolamine synthesis, membrane biology, and neurodegenerative diseases.

[Learn More](#)

Parameter	Specification
Product Name	PCYT2 Rabbit pAb
Remarks/Alias	ET, SPG82, PCYT2
Species	Human
Gene ID (Human)	5833
Gene ID	5833
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-389 of human PCYT2 (NP_002852.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99447
Protein Weight	44kDa
Shipping	Ice bag