

Anti-Pd-1/Cd279 Polyclonal Antibody For Wb, If-P, Elisa - Q02242

Item Number: CM0022636



Introduction

Rabbit polyclonal antibody against mouse PD-1/CD279, an inhibitory immune checkpoint receptor, suitable for Western blot, immunofluorescence, and ELISA. Ideal for studying T-cell tolerance and tumor immunology.

[Learn More](#)

Parameter	Value
Product Name	PD-1/CD279 Rabbit pAb
Remarks/Alias	PD-1; Pdc1; Ly101; PD-1/CD279
Species	Mouse
Gene ID	18566 (Mouse)
Immunogen	Recombinant protein of mouse PD-1/CD279
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse
SWISS	Q02242
Protein Weight	32kDa
Shipping	Ice bag