

Anti-Pd-1/Cd279 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q15116

Item Number: CM0018677



Introduction

Rabbit polyclonal antibody against human PD-1/CD279 (SWISS Q15116), validated for WB, IF-P, ELISA. Cross-reacts with mouse and rat. Ideal for immuno-oncology research and IVD development.

[Learn More](#)

Parameter	Value
Product Name	PD-1/CD279 Rabbit pAb
Remarks/Alias	PD1; PD-1; CD279; SLEB2; hPD-1; hPD-I; hSLE1; PD-1/CD279
Species	Human, Mouse, Rat
Gene ID	5133
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 205-259 of human PD-1/CD279 (NP_005009.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15116
Protein Weight	32kDa
Shipping	Ice bag