

Anti-Pd-L1/Cd274 Monoclonal Antibody For Fc - Q9Nzq7

Item Number: CM0005952



Introduction

-conjugated rabbit anti-human PD-L1/CD274 monoclonal antibody validated for flow cytometry. Reacts with human. Explore PD-L1/PD-1 pathway in immune checkpoint research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human PD-L1/CD274 mAb
Remarks/Alias	B7-H; B7H1; PDL1; PD-L1; hPD-L1; PDCD1L1; PDCD1LG1
Species	Human
GeneID (human)	29126
GeneID	29126
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human PD-L1/CD274 (Q9NZQ7
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9NZQ7
Protein Weight	33kDa
Shipping	Ice bag