

Anti-Pd-L1/Cd274 Monoclonal Antibody For Wb, Elisa - Q9Nzq7

Item Number: CM0002284



Introduction

Recombinant rabbit monoclonal antibody against human PD-L1/CD274 (Q9NZQ7), suitable for Western blot and ELISA. Essential for research on immune checkpoint regulation, T-cell modulation, and tumor immunoescape mechanisms.

[Learn More](#)

Parameter	Value
Product Name	PD-L1/CD274 Rabbit mAb
Remarks/Alias	B7-H; B7H1; PDL1; PD-L1; hPD-L1; PDCD1L1; PDCD1LG1; PD-L1/CD274
Species	Human
GeneID (Human)	29126
GeneID	29126
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human PD-L1/CD274 (Q9NZQ7)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NZQ7
Protein Weight	33kDa
Shipping	Ice bag