

Anti-Pd-L1/Cd274 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa, Inhibition - Q9Nzq7

Item Number: CM0011158



Introduction

PD-L1/CD274 rabbit monoclonal antibody raised against recombinant protein, for WB, IF/ICC, ELISA, and inhibition assays. Specific for human PD-L1, validated for immune checkpoint research and cancer immunotherapy studies.

[Learn More](#)

Parameter	Value
Product Name	PD-L1/CD274 Rabbit mAb
Remarks/Alias	B7-H; B7H1; PDL1; PD-L1; hPD-L1; PDCD1L1; PDCD1LG1
Species	Human
Gene ID (Human)	29126
Gene ID	29126
Immunogen	Recombinant protein of human PD-L1/CD274
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA, Inhibition
Cross Reactivity	Human
SWISS	Q9NZQ7
Protein Weight	33kDa
Shipping	Ice bag