

Anti-Pd-L2/Cd273/B7-Dc Rabbit Polyclonal Antibody For Wb, Elisa - Q9Wul5

Item Number: CM0022870



Introduction

Rabbit polyclonal antibody targeting mouse PD-L2 (CD273/B7-DC), suitable for Western blot and ELISA. Cross-reacts with rat. Useful for research in immune tolerance and checkpoint pathways.

[Learn More](#)

Parameter	Value
Product Name	PD-L2/CD273/B7-DC Rabbit pAb
Remarks/Alias	Btdc; B7-DC; PD-L2; F730015O22Rik; PD-L2/CD273/B7-DC
Species	Mouse
GeneID (human)	58205
GeneID	58205
Immunogen	Recombinant protein of mouse PD-L2/CD273/B7-DC.
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9WUL5
Protein Weight	28kDa
Shipping	Ice bag