

Anti-Pdc Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P20941

Item Number: CM0025072

Introduction

Anti-PDC Rabbit Polyclonal Antibody for detection of phosducin (28 kDa) in WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Validated for visual phototransduction and photoreceptor metabolism studies. Research use.

[Learn More](#)

Parameter	Value
Product Name	PDC Rabbit pAb
Remarks/Alias	PHD; MEKA; PhLP; PhLOP; PDC
Species	Human
Gene ID (Human)	5132
Gene ID	5132
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-246 of human PDC (NP_002588.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20941
Protein Weight	28kDa
Shipping	Ice bag