


Anti-Pdcd10 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Bul8

Item Number: CM0019634



Introduction

Rabbit polyclonal antibody against human PDCD10 (CCM3/TFAR15) for Western blot, IHC-P, IF/ICC, and ELISA. Part of STRIPAK complexes, involved in apoptosis, angiogenesis, and cell migration. Cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	PDCD10 Rabbit pAb
Remarks/Alias	CCM3, TFAR15, PDCD10
Species	Human
GeneID (Human)	11235
GeneID	11235
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-212 of human PDCD10 (NP_009148.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BUL8
Protein Weight	25kDa
Shipping	Ice bag