


## Anti-Pdcd4 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q53EI6

Item Number: CM0009443



### Introduction

PDCD4 Rabbit mAb for WB, IHC-P, ELISA applications. Recognizes human, mouse, rat PDCD4 (52 kDa), a translation inhibitor and tumor suppressor. Ideal for cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                                   |
|------------------|-------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PDCD4 Rabbit mAb                                                                                                        |
| Remarks/Alias    | H731; PDCD4                                                                                                             |
| Species          | Human                                                                                                                   |
| GeneID (human)   | 27250                                                                                                                   |
| GeneID           | 27250                                                                                                                   |
| Immunogen        | Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 370-469 of human PDCD4 (Q53EL6). |
| Source           | Rabbit                                                                                                                  |
| Category         | Monoclonal Antibodies                                                                                                   |
| Application      | WB, IHC-P, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                       |
| SWISS            | Q53EL6                                                                                                                  |
| Protein Weight   | 52kDa                                                                                                                   |
| Shipping         | Ice bag                                                                                                                 |